


# PRODUCT DATA SHEET

Revision 1	Cell Scraper	
Validfor Item-No.:	5050000/5052400/5052000/5054400	

1.	Description / Specification	
1.1	Description	Cell Scraper, sterile 5050000: Cell Scraper, 18 cm, Blue, Blade and Handle Level, Sterile 5052400: Cell Scraper, 23 cm, yellow, blade mounted perpendicular to handle, sterile 5052000: Cell Scraper, 25 cm, green, blade and handle level, sterile 5054400: Cell Scraper, 30 cm, yellow, blade and handle placed vertically, sterile
1.2	Dimensions	see Customer Drawing
1.3	Volume	NA
1.4	Material / Resin	<u>Handle:</u> PP (Polypropylene) <u>Blade:</u> TPU(Thermoplastic Polyurethane) The materials for manufacturing are <i>USP Class VI</i>
1.5	Colour	<u>Handle:</u> 5050000 Blue, 5052400/5054400 Yellow, 5052000 Green <u>Blade:</u> White
1.6	Sterilization	sterile
1.7	Quality Control	- <u>Raw Material-Control:</u> physical testing - <u>Product-Control:</u> testing of attributive and variable characteristics in accordance with the valid specification
1.8	Intended Use	General laboratory product for bacteriology to be used by qualified personnel in a laboratory environment.
1.9	Other Information	- For single use only

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic
2.2	Temperature range	For application: -20°C to +60 °C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Shelf life	3 years
2.6	Other Information	-

3.	Packaging	
3.1	Pieces / Bag	1
3.2	Pieces / Box	100
3.3	Lot-No.	DDD YY XXX (day, year, product No.)
3.4	Other Information	Certificate of Quality to download

4.	Other Information
4.1	Laboratory use only

Prior Issue	Drawn	Approved	Released	<b>CONFIDENTIAL:</b> Information contained in this document or drawing is confidential and proprietary to Zhejiang Saining Biotechnology Co., Ltd. This document may not be reproduced for any reason without written permission from Zhejiang Saining Biotechnology Co., Ltd. All rights of design, invention, and copyright are reserved.
Revision 1	Date 1 June 2024	Date 1 June 2024	Date 1 June 2024	
Date 01.06.2024	Name Chen Anjing	Name Zhang Heng	Name Wang Kedong	

**DISCLAIMER:** The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.

